

Name: Jannatul Ferdous Umama


ID: 20-42616-1

Sec: R



Experiment - 6: Determination of Ferrous Ion (Fe^{2+}) in a supplied of Iron salt by standard Potassium Dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$) Solution.

Experimental data:

Method: Oxidation - Reduction titration



Redox Half reactions:



Weight taken (in gm) = 0.66 g

The strength of $\text{K}_2\text{Cr}_2\text{O}_7$ Solution

$$= \frac{\text{weight taken (in gm)} \times 0.1}{0.49} \text{ (N)}$$

$$= \frac{0.66 \times 0.1}{0.49} \text{ (N)}$$

$$= 0.1346 \text{ (N)}$$

$$= 0.13 \text{ (N)}$$

Table : Determination of the amount of iron in Mohr's salt solution using standard $K_2Cr_2O_7$ Solution.

No of reading	Vol. of Mohr's Salt Solution	Vol of $K_2Cr_2O_7$ (burette reading) (in mL)			Mean (in mL) (v)
		Initial	Final	Difference	
1	10	0.00	5.00	5.00	5.03
2	10	5.00	9.90	4.90	
3	10	9.90	15.00	5.1	
4	10	15.00	20.10	5.1	

Calculations: $1\text{ mL } 1N K_2Cr_2O_7 \equiv 0.05584\text{ gm of } Fe^{2+}$

Amount of iron in 10mL of iron salt solution

$$0.05584 \times V \times S \text{ (gm)}$$

$$= 0.05584 \times 5.03 \times 0.13$$

$$= 0.0365\text{ gm}$$

$$= 0.04\text{ gm}$$

Amount of iron in 500mL of iron salt solution

$$0.05584 \times V \times S \times 50\text{ gm}$$

$$= 0.05584 \times 5.03 \times 0.13 \times 50$$

$$= 1.8256\text{ gm} = 1.83\text{ gm}$$

Observe value of Fe^{2+} (in 500mL solution)

$$= 1.83\text{ gm}$$

$$\text{Known value of Fe}^{2+} (\text{in } 500\text{ml solution}) \\ = \frac{55.84 \times 8.55}{392.14} (\text{gm})$$

$$= 1.2175 (\text{gm})$$

$$= 1.22 \text{ gm.}$$

Result : The amount of Ferrous ions in 500ml of iron salt solution is 0.

Percentage of error :

$$\frac{\text{Known value} - \text{Observed value}}{\text{Known value}} \times 100 \\ = \frac{1.22 - 1.83}{1.22} \times 100 \\ = -0.28\%$$