

Exp 1. Standardization of Sodium Hydroxide (NaOH) Solution with Standard oxalic Acid ($C_2H_2O_4 \cdot 2H_2O$) Solution.

Name : Jannatul Ferdows Umama , ID: 20-42616-1
Section : R.

Method : Acid based titration

Experimental Data :

$$\begin{aligned} &\text{Strength of oxalic acid solution} \\ &= \frac{\text{Weight taken (in gm)} \times 0.1}{0.63} \text{ (N)} \\ &= \frac{0.63 \times 0.1}{0.63} \text{ N.} \\ &= 0.1 \text{ N} \end{aligned}$$

Table : Standardization of supplied NaOH solution against standard oxalic acid solution by acid-base titration

No of reading	Vol. of NaOH (in ml)	Vol. of Oxalic acid (burette reading) (in ml)			Mean (in ml)
		Initial	Final	Difference	
1	10	0.00	8.50	8.50	8.50 +
2	10	8.50	17.10	8.60	8.60 +
3	10	17.10	25.60	8.50	8.50
					3
					= 8.533

③ Calculations :

The strength of supplied NaOH Solution :

$$V_{\text{NaOH}} \times N_{\text{NaOH}} = V_{\text{oxalic acid}} \times N_{\text{oxalic acid}}$$

$$N_{\text{NaOH}} = \frac{8.53 \times 0.1}{10} = 0.0853 \text{ N}$$

$$= 0.09 \text{ (N)}$$

Results : The strength of Supplied NaOH solution is 0.09 (N).